

Standard operating procedure of Nanodrop

1.Objective: This instruction explains how use and keep the Nanodrop

2.Scope: using for documenting the concentration of nucleic acids or proteins

3.Responsibility: Laboratory assistance is responsible for the accuracy of the apparatus and should be reported to him/her in case of problems.

4. Materials and equipment

5.Guide lines for use:

Turn on the Nanodrop and its computer.

Run the ND-1000 software on the Based on your sample type, select Nucleic acid or protein options.

Wash the sensor of the apparatus with 2 μ L deionized water (DW) and whatman filter paper

Again, add 2 μ L deionized water on the center of sensor, lower down the arm of Nanodrop and press the “blank” key in the software.

Again remove the DW and add 2 μ L of your source solution that used for extraction procedure - ∞ on the sensor and press the Blank key.

After removing the source solution, add 2 μ L of your sample and lower down the arm and - ∞ press the “Measure” key.

Document the results and after washing the sensor with DW, exit from the software and turn - ∞ off the apparatus.

6.Warning and safety precaution:

If any abnormal sound or vibrate are felt, press the “Stop” bottom immediately

7.Maintenance:

The door of the apparatus should be always closed.

Document ID: CoreLab. 15

Document name	Standard operating procedure of Nanodrop
Date	99/4/1
Full file name	Standard operating procedure of Nanodrop
Document description	These instructions describe the usage, maintenance and calibration of Nanodrop which is kept in Core laboratory.
Producer	Core Lab
