

# cDNA synthesis (RT-PCR)

#### Workflow

RNA Extraction ———— cDNA synthesis (RT-PCR) ———— Real-Time PCR

cDNA synthesis, also known as reverse transcription, generates complementary DNA (cDNA) from an RNA template.

### **Specimen**

RNA template

### **Materials**

Random hexamer primer

Reverse transcriptase

RNase inhibitor

dNTP

microfuge tube

Ultrapure DNase/RNase

Free Distilled water

## **Equipment and supplies**

Refrigerated Microcentrifuge

Centrifuge tubes

Pipette tips

Nanodrop

Vortex Mixer

Thermocycler

water bath



### Safety (Warning and Biohazard consideration):

Microcentrifuge tubes and Pipette tips must be autoclaved before use. A major source of RNase contamination is form the hands of the researcher. Gloves should be worn all stages during the preparation of materials and solutions used for the isolation and analysis of RNA, and during all manipulations involving RNA.

#### **Procedure:**

**Step1**. Prepare the following mixture in a microtube.

Reagent	volume	
Oligo dT Primer (50 μM)	1 μl	
or Random 6 mers (50 μM)	or 1 μl (0.4 - 2 μl)	
dNTP Mixture (10 mM each)	1 μl	
Template RNA	Total RNA: < 5 μg	
RNase Free dH2O PolyA+ RNA: < 1 μg		
	x μl	
total	10 μl	

Step2. Incubate for 5 min at 65°C, then cool immediately on ice

Step3. Prepare the reaction mixture in a total volume of 20 µl.

Reagent	volume	
Template RNA Primer Mixture (from step 2)	10 μl	
RNase Inhibitor (40 U/μl)	0.5 μl (20 U)	
5X PrimeScript Buffer	4 μl	
PrimeScript RTase (200 U/μl)	1.0 µl (200 U)	
RNase Free dH2O	x μl	
Total	20 μl	

**Step4**. Mix by pipetting up and down.

Add Buffer /Enzyme mix to the first tube, mix gently and transfer in a 0.2 mL microfuge tube.

**Step5**. Incubate the reaction mixture using the following conditions.

30°C 10 min (required when using Random 6 mers)

42°C (50°C) 30 - 60 min 95°C 5 min



Step 6. Analyze the PCR products by agarose gel electrophoresis

Step 7. Store cDNA at -20°C

### **References:**

 $\underline{https://www.neb.com/protocols/1/01/01/first-strand-cdna-synthesis-e6300}$ 

https://www.takara-bio.com

Core 1	Lab.105	به سند:	شناسناه
	Luvitus		·

 نام سند
 پروتکل استاندارد سنتز CDNA

 تاریخ صدور
 ۱۳۹۹/7/7

 نام کامل فایل
 پروتکل سنتز CDNA

 شرح سند
 روش سنتز CDNA را شرح می دهد.

 تهیه کننده
 شهربانو نادری